# Facebow®

## Patient Instructions

Series 105

### Facebow Patient Instructions

Your doctor has prescribed the use of the enclosed facebow for your orthodontic treatment. It is very important for you to follow your doctor's directions on the use and care of this device.

Improper or careless use can be dangerous and could cause serious injury.

### Important Safety Instructions

1. Carefully follow your doctor's wearing instructions.
2. Never place or remove your facebow unless the headgear straps have been disconnected.
3. Do not wear your facebow while participating in sports or running.
4. Do not allow anyone to grab or pull on your facebow or headgear.
5. Never lift the facebow over your face.
6. Do not adjust the facebow without first disconnecting the straps.
7. Call your doctor if you ever have any questions on the proper use of your facebow.

### CAUTIONS

Improper or careless use can be dangerous and could cause serious injury. Failure to specifically follow your doctor's instructions could result in serious injury, including but not limited to, damage to your eyes and/or blindness.

This product is manufactured with nickel which may cause an allergic reaction in some people. Discontinue use immediately and call your orthodontist if any reaction occurs.

**NOTE:** This is a single-use device. Any reuse of this product will run the risk of cross-contamination to patients.

### NOTICE:

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

Order online at [tportho.com](http://tportho.com).

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