INSTRUCTIONS

Free-N-Easy[®] MouthGuard

Series 100

 Made by
TP Orthodontics, Inc.
100 Center Plaza
La Porte, Indiana 46350-9672 USA

TPOrtho.com

 Auth. EU Rep.
TP Orthodontics Europe
Fountain Ct.
12 Bruntcliffe Way
Morley, Leeds LS27 0JG UK

Free-N-Easy is a registered trademark of TP Orthodontics, Inc.

© 2006 TP Orthodontics, Inc.

996-520, Rev. 6/14

ENGLISH

Free to talk, easy to breathe

You're free to talk, and it's easy to breathe when you wear this new, improved hinge-type mouth guard.

The Free-N-Easy mouth guard protects both the upper and lower teeth and lips—even when your mouth is open.

Molded with smooth surfaces, from soft plastic, the Free-N-Easy mouth guard is simple to clean and provides maximum comfort.

- The sidewalls do not extend to end of hinges, so gag reflex and "full feeling" are reduced
- Integral hinge permits talking and facilitates breathing while playing
- Interlocking serrations on upper and lower distribute forces evenly to both jaws
- Dimensional to fit both children and adults
- Will not interfere with orthodontic appliances
- Molded from soft, durable plastic
- All edges rounded for maximum comfort
- Available with safety straps for contact sports

Free-N-Easy Mouth Guard Care Instruction

Wash before using and after each use with warm water (do not boil). Rinse and dry. It is important to follow these instructions to reduce the possibility of oral disease or infection.

CAUTION: The Free-N-Easy mouth guard has been designed to reduce the incidence of injuries to the mouth in connection with certain contact sports. However, no mouth guard is totally effective in eliminating the incidence of injuries.

Remember—sound, healthy teeth are more resistant to trauma. It is suggested that every player have a dental examination.

NOTICE:

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

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