

LIGATURE COLOR SELECTOR*



TP Orthodontics, Inc.
tportho.com



Pale Pink



Lavender



Mint



Tooth-Colored



Clear



Pink



Purple



Green



Blue



Black



Red



Yellow



Orange



Light Blue



White



Silver-Gray



Evergreen



Navy Blue



Burgundy



Teal

Mini Stix Ligature Ties

Color	Super Slick®	Uncoated
Pale Pink	382-916	383-916
Lavender	382-917	383-917
Mint	382-918	383-918
Tooth-Colored	382-919	383-919
Clear	382-921	383-921
Pink	382-922	383-922
Purple	382-923	383-923
Green	382-924	383-924
Blue	382-925	383-925
Black	382-927	383-927
Red	382-928	383-928
Yellow	382-929	383-929
Orange	382-930	383-930
Light Blue	382-931	383-931
White	382-932	383-932
Silver-Gray	382-934	383-934
Evergreen	382-937	383-937
Navy Blue	382-938	383-938
Burgundy	382-939	383-939
Teal	382-940	383-940
Assorted Colors - 5 of each color		383-941

Super Slick® Mini Stix are virtually frictionless elastomeric ligature ties. They feature a highly lubricated Metafasix® technology coating that, when moistened, provides a slippery surface that significantly reduces friction at the point of archwire/ligature contact.

Uncoated Mini Stix elastomers are designed for single patient use to reduce the risk of cross-contamination. The elastomeric material used for manufacturing is highly resilient, colorfast and stain-resistant.

Mini Stix Ligature Ties are available in packages of 40 stix. The assorted colors have 100 stix, 5 of each color. Mini Stix are latex-free.

Visit **tportho.com** to order ligatures or call **800-348-8856**.



TP Orthodontics, Inc.
tportho.com

*Colors may vary.

Super Slick is a registered trademark of TP Orthodontics, Inc. and manufactured under US Patent 6,203,317; Britain 1174107; France 1174107; Germany 60113016.2; Italy 1174107. Metafasix is a registered trademark of TP Orthodontics, Inc. and manufactured under US Patent 6,203,317.

© 2008 TP Orthodontics, Inc.

996-285ENGU, Rev. 0