


Nu-Edge[®] lite

Cobalt Chromium Pre-Coated Brackets

- Aesthetic polymer mesh base.**
- Biocompatible for nickel sensitive patients.**
- Pre-coated for convenience.**




The Nu-Edge lite pre-coated bracket combines the performance of cobalt chromium with the aesthetics of a clear polymer mesh base. The result is superior performance with enhanced aesthetics.



Exceptional aesthetics increases patient satisfaction with small design and clear base for a minimal visual footprint.

- Strong cobalt chromium; 100% nickel free
- Clear base allows light to penetrate for a full light cure, resulting in improved bond strength
- Increases ease of application and reduces chair time
- Simply remove bracket from the package with bonding forceps, place and cure
- Available in MBT prescription



Save Time with Pre-Applied Adhesive

Nu-Edge lite Brackets feature Readi-Base[®], with eXact[®] pre-applied adhesive for convenience and increased efficiency.

- Ready to place right out of the package
- Reduced risk of contamination
- Brackets set easily without sliding
- Minimal flash clean-up



TP Orthodontics, Inc.
tportho.com

Nu-Edge® lite Bracket

Cobalt Chromium with Pre-Applied Adhesive

MBT* McLaughlin, Bennett, Trevisi System

Patient Kits

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Upper & Lower 5-5	Cuspid/Bicuspid	285-001RBX	285-004RBX

Kits contain 20 all metal brackets.

Packs of 5 — Maxillary Brackets

	Slot Size	Left/ Right	Bondable
Central Torque +17° Tip +5°	0.46mm (.018")	L	285-010RBX
		R	285-011RBX
Lateral Torque +10° Tip +7°	0.46mm (.018")	L	285-014RBX
		R	285-015RBX
Cuspid w/ Distal Hook Torque 0° Tip +7°	0.46mm (.018")	L	285-022RBX
		R	285-023RBX
Universal Bicuspid w/ Hook Torque -7° Tip 0°	0.46mm (.018")	L	285-030RBX
		R	285-031RBX
	0.56mm (.022")	L	285-012RBX
		R	285-013RBX
	0.56mm (.022")	L	285-016RBX
		R	285-017RBX
	0.56mm (.022")	L	285-024RBX
		R	285-025RBX
	0.56mm (.022")	L	285-032RBX
		R	285-033RBX


Package of 5.

Packs of 5 — Mandibular Brackets

	Slot Size	Left/ Right	Bondable
Anteriors Torque -6° Tip 0°	0.46mm (.018")	L/R	285-054RBX
		L/R	285-055RBX
Cuspid w/ Distal Hook Torque 0° Tip +3°	0.46mm (.018")	L	285-060RBX
		R	285-061RBX
1st Bicuspid w/ Hook Torque -12° Tip +2°	0.46mm (.018")	L	285-068RBX
		R	285-069RBX
2nd Bicuspid w/ Hook Torque -17° Tip +2°	0.46mm (.018")	L	285-076RBX
		R	285-077RBX
	0.56mm (.022")	L	285-062RBX
		R	285-063RBX
	0.56mm (.022")	L	285-070RBX
		R	285-071RBX
	0.56mm (.022")	L	285-078RBX
		R	285-079RBX

Package of 5.

MBT System



996-478ENG, Rev. 1

Readi-Base is a registered trademark of TP Orthodontics, Inc. and manufactured under US Patents 6,685,468, 6,834,761, 7,131,836, 7,780,442, RE45,743; Britain 1562511; China ZL 03824889.1, 03822189.6; Europe 1545380; France 1562511; Germany 60332263.8-08; Italy 1562511; Japan 4503440, 4414334; Mexico 274623, 265725; Spain 1562511, ES2362089T3. Nu-Edge and eXact are registered trademarks of TP Orthodontics, Inc. © 2017 TP Orthodontics, Inc. All rights reserved.

