

Micro6™ Molar Tube

with Zip-Release

Unique zip-release cap converts only when you want it to.

TPO® announces a new Single Convertible Molar Tube — Micro6 with Zip-Release. The new upper bondable tube has a smaller profile height with fully functional tie wings. The convertible cap is a newly designed zip-release cap for easier conversion.



Features

- Low profile height
- Fully functional tie wings for easy ligation
- Rounded ball-end hooks
- Perforated cap for easy conversion
- Larger fluted tube opening
- Cobalt chromium for strength and biocompatibility
- Anatomically contoured base

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset
0.46 x 0.64mm (.018" x .025")	4.0	-10°	0°	4°
0.56 x 0.71mm (.022" x .028")	4.0	-10°	0°	4°

A standard cap removing plier can be used to convert the Micro6 cap.

The Micro6 tube is an addition to the convertible tube product line and for doctors that prefer to treat with a -10 torque.

Find out more at tportho.com

