

# **Surgical Supplies for the Orthodontic Patient**

From TP Orthodontics



TP Orthodontics, Inc.



**TP Orthodontics,  
a leading  
manufacturer  
of orthodontic  
supplies, offers a  
complete line of  
traction hooks  
and accessories  
for orthodontic  
patients.**

## Bondable Traction Hooks

Designed to encourage and direct the eruption of malpositioned or impacted teeth when the crown is surgically exposed, TPO's traction hooks with a chain can be used when the tooth is deeply embedded in the gingival tissue.



T P O R T H O . C O M

# Round and Rumely Traction Hooks



**Round Traction Hook**




**Attach the traction hook to the orthodontic archwire with Zing String elastomeric thread.**



**Round Traction Hook with Stiff Chain**

The stiff ligating chain is used when a tooth is deeply embedded in gingival tissue.



**Stiff ligating chain can be cut to length.**



**Rumely Traction Hook with Flexible Chain**



**Links are sized for easy engagement of elastomeric thread or wire.**

# Elastomeric Thread

## Slip Free Zing String®

- Aesthetic material
- Special cohesive surface helps keep knots tied
- Use for attaching to traction hooks for orthodontic traction




## Crimpable Hooks



Inner surface features an abrasive coating that results in force loads equal to welded hooks.

## Ligature Wires

Preformed; length is ample for use with all tying pliers.



## Bonding Preparation

### TPO® Blue Etchant Gel

- Can be used with any bonding system
- Pre-loaded in a syringe
- Includes supply of disposable tips
- Curved tip simplifies application on high cuspids, bicuspids, and molars
- Etches enamel in 30 seconds
- Phosphoric acid content 37%



### eXact®

## Orthodontic Adhesive

- Medium-high viscosity light cure adhesive
- Use for metal, plastic, or ceramic attachments
- Prevents attachment drift, produces little flash, and is extremely easy to clean up



# TP Orthodontics

## Surgical Components

### Parts List

- 224-011** Round Traction Hook (10ct)
- 224-012** Round Traction Hook with Stiff Chain (10ct)
- 224-014** Rumely Traction Hook with Flexible Chain (10ct)
- 100-247** Slip-free Zing String; Light (.018") diameter (25ft)
- 100-245** Slip-free Zing String, Medium (.020") diameter (25ft)
- 100-246** Slip-free Zing String, Heavy (.025") diameter (25ft)
- 382-001** Stainless Steel .025mm ligature wire (1000ct)
- 226-010** Crimpable Hooks (10ct)
- 100-172** Crimpable Hook Placement Pliers
- 100-132** Straight Tip Hemostat
- 151-270** eXact Light Cure Adhesive Kit
- 151-271\*** eXact Light Cure Adhesive in Syringe, 4g (2ct)
- 150-713\*** TPO Blue Etchant Gel in Syringe (5g, 10 tips)
- 150-714** Etchant Gel Disposable Dispensing Tips (25ct)
- 151-272\*** Light Cure Sealant (5cc)
- 150-704\*** Disposable Brushes (50ct)

*\*Components of eXact Light Cure Adhesive Kit, part 151-270.*

**Orthodontists around  
the world trust TPO®  
to provide efficient  
and predictable  
treatment options.**



TP Orthodontics, Inc.

100 Center Plaza  
La Porte, Indiana 46350-9672  
United States  
800.348.8856  
info@tportho.com



T P O R T H O . C O M